


DEX™

A WIRELESS WORLD


COM-DEX

UK
DE
FR
ES
DK
IT
JP
US
CA


WIDEX®
HIGH DEFINITION HEARING

4

5

UK	USER INSTRUCTIONS.....	2
DE	BEDIENUNGSANLEITUNG.....	16
FR	MODE D'EMPLOI	30
ES	INSTRUCCIONES DE USO	44
DK	BRUGSANVISNING	58
IT	LIBRETTO D'ISTRUZIONI.....	72
JP	取扱説明書.....	86
US/CA	USER INSTRUCTIONS.....	100

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Intended use

The device is intended to improve the ability of hearing impaired people to conduct mobile phone conversations while using their hearing aids, to improve the usability and ease of control of such hearing aids by equipping them with a Smartphone-based control method.

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2.1 USB cable

CHARGING – *See illustration 3*

- 3.1 Connect USB cable between the device and the power supply. We recommend that the two neck loop parts are not joined while charging the device.
- 3.2 Red light: charging, green light: fully charged.

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How to pair the device with your phone

To pair your device with your phone for the first time, follow these steps:

1. Make sure that the device and phone are within 1 metre range of each other.
2. Connect the neck loop. The device will blink red and blue.
3. Activate Bluetooth on your phone.
4. Start the pairing process according to the instructions in your phone's user guide.
5. Your phone will show a list of Bluetooth devices.
Select COM-DEX from the list.

Once you have paired the device, it will remain so even when you disconnect the neck loop.

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Note: Answering and ending calls is also possible on the mobile phone. Note that, on some mobile phones, you may need to select the COM-DEX from a list of devices. Otherwise the sound will be transmitted through the phone's speakers.

Outgoing calls – *See illustration 5*

- 5.4 Dial from phone.

Note: You can also use the device to stream sound from any sound source to your hearing aids.

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You can download the COM-DEX app from the App Store or Google Play. This app allows you to use your phone as a remote control.

MAINTENANCE

- Clean with a damp cloth. Do not use chemicals.
- Do not expose to extreme temperatures or high humidity.
- Do not immerse in liquids.
- When the device is not in use keep it in a dry location out of reach of children and pets.
- Never try to open or repair the device yourself. This may only be done by authorised personnel.

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TROUBLESHOOTING

Problem	Possible cause	Solution
The device does not work	a) Device battery low b) The mobile phone is set to silent mode c) Hearing aid battery low d) Out of range	a) Charge device battery b) Set to normal mode c) Change hearing aid's battery d) Bluetooth has a range up to 10 m.
No sound	a) Not connected via Bluetooth	a) Repeat the pairing process

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SYMBOLS

Symbols commonly used by Widex A/S in medical device labelling (labels/IFU/etc.)

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Symbol	Title/Description
	<p>Catalogue number The product's catalogue (item) number.</p>
	<p>Serial number The product's serial number.*</p>
	<p>Keep away from sunlight The product must be protected from light sources and/or The product must be kept away from heat</p>
	<p>Keep dry The product must be protected from moisture and/or The product must be kept away from rain</p>
	<p>Lower limit of temperature The lowest temperature to which the product can be safely exposed.</p>

Symbol	Title/Description
	<p>Upper limit of temperature</p> <p>The highest temperature to which the product can be safely exposed.</p>
	<p>Temperature limits</p> <p>The highest and lowest temperatures to which the product can be safely exposed.</p>
	<p>Consult instructions for use</p> <p>The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.</p>
	<p>Caution/Warning</p> <p>Text marked with a caution/warning symbol must be read before using the product.</p>

Symbol	Title/Description
	<p>WEEE mark “Not for general waste”</p> <p>When the product is to be discarded, it must be sent to a designated collection point for recycling and recovery.</p>
	<p>CE mark</p> <p>The product is in conformity with the requirements set out in European CE marking directives.</p>
	<p>Alert</p> <p>The product is identified by R&TTE Directive 1999/5/EC as an equipment Class 2 product with some restrictions on use in some CE member states.</p>
	<p>RCM mark</p> <p>The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand market.</p>

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DESCRIPTION OF DEVICE – *See illustration 1*

- 1.1 Push button
- 1.2 Room Off button
- 1.3 LED
- 1.4 Microphone
- 1.5 Neck loop
- 1.6 Micro USB

Note: Read these instructions carefully before you start using the device.

Intended use:

This device is intended to improve your ability to use Bluetooth enabled mobile phone connections with your hearing aid(s). The aid(s) can further be manipulated by the controls on the mobile phone.

Note: The device uses Bluetooth standard which is common in newer mobile phones.

ACCESSORIES – *See illustration 2*

2.1 USB cable

CHARGING – *See illustration 3*

- 3.1 Connect USB cable between the device and any power supply. We recommend that the two neck loop parts are not joined while charging the device.
- 3.2 Red light: charging, green light: fully charged.

Battery life: 8 hours streaming, 8 days standby.

HOW TO USE – *See illustration 4*

- 4.1 Join the two neck loop parts. The neck loop functions as an on/off switch.
- 4.2 Pair the device to your mobile phone using Bluetooth.

How to pair the device with your mobile phone

To pair your device with your mobile phone for the first time, follow these steps:

1. Make sure that the device and mobile phone are within 3 ft range of each other.
2. Connect the neck loop. The device will blink red and blue.
3. Activate Bluetooth on your mobile phone.
4. Start the pairing process according to the instructions in your mobile phone's user guide.
5. Your mobile phone will show a list of Bluetooth devices. Select COM-DEX from the list.

Once you have paired the device, it will remain paired even when you disconnect the neck loop.

Incoming calls – *See illustration 5*

- 5.1 Answer call: Press push button.
- 5.2 Speak into microphone while device is hanging around your neck.
- 5.3 End call: Press push button.

Note: Answering and ending calls is also possible on the mobile phone. Note that, on some mobile phones, you may need to select the COM-DEX from a list of devices. Otherwise the sound will be transmitted through the mobile phone's speakers.

Outgoing calls – *See illustration 5*

- 5.4 Dial from mobile phone.

Note: You can also use the device to stream sound from any Bluetooth enabled sound source to your hearing aids.

Sound adjustment during use

Volume is adjusted on the mobile phone.

To shut out surrounding sounds: press Room Off button – See *illustration 1*.

You can download the COM-DEX app from the App Store or Google Play. This app allows you to use your mobile phone as a remote control.

MAINTENANCE

- Clean with a damp cloth. Do not use chemicals.
- Do not expose to extreme temperatures or high humidity.
- Do not immerse in liquids.
- When the device is not in use keep it in a dry location out of reach of children and pets.
- Never try to open or repair the device yourself. This may only be done by authorized personnel.

INFORMATION ABOUT EMISSION

Emission levels from Widex devices are several times lower than the levels transmitted by mobile phones and other communications equipment.

For further information, especially regarding patients with medical implants, download the document: “Information on WidexLink™” at www.widex.com/widexlink

TROUBLESHOOTING

Problem	Possible cause	Solution
The device does not work	a) Device battery low b) The mobile phone is set to silent mode c) Hearing aid battery low d) Out of range	a) Charge device battery b) Set to normal mode c) Change hearing aid battery d) Bluetooth has a range up to 30 ft.
No sound	a) Not connected via Bluetooth	a) Repeat the pairing process

If the problems persist, contact your hearing care professional.

WARNINGS

WARNING

- Do not use the device in mines or other areas with explosive gases.
- When the device is connected via micro USB cable to external mains-operated equipment such as a charger, PC, audio source or similar, this equipment must comply with IEC 60065, IEC 60950-1, IEC 60601-1 or equivalent safety standards.

REGULATORY INFORMATION

The following Table summarizes the technical details of the Bluetooth and WidexLink technology as it is implemented in the COM-DEX.

Unit	COM-DEX		Hearing aid
Type	W-link	Bluetooth*	W-link
Antenna type	Inductive Loop antenna	Embedded ceramic antenna	Inductive antenna
Antenna dimensions	L: 70 cm	NA	Ø1.8 mm L – 4.85mm
Modulation	FSK	FHSS/GFSK, $\pi/4$ DPSK, 8DPSK	FSK

Magnetic Field Strength	-26dB μ A/m at 10m distance	NA	-54dB μ A/m at 10m distance
Output power	1.2 nW EIRP** -59 dBm	2.5 mW EIRP +4dBm	29pW EIRP** -75 dBm
Range	<1m remote unit to hearing aid	<10m between COM-DEX and Bluetooth device	<1m remote unit to hearing aid <30cm between hearing aids
Center frequency	10.6MHz	2.45GHz ISM	10.6MHz
Channels	Single channel radio	80 channels	Single channel radio

Bandwidth	660kHz (-15dB)	1MHz	660kHz (-15dB)
Duty Cycle ***	Up to 100 %	< 100 %	< 5 %
Data rate	212kbit/s raw channel capac- ity	3 Mbit/s	212kbit/s raw channel capac- ity
Data flow	Simplex or semiduplex	Time division duplex (TDD)	Simplex or semiduplex
Protocol	Random access - no collision avoidance	Packet based protocol, time divided	Random access - no collision avoidance

** Bluetooth specification v 4.1 + EDR/LE published by the Bluetooth Special Interest Group (SIG)*

*** EIRP = Equivalent isotropically radiated power*

**** Duty Cycle averaged over 1 hour of operation*

Cables and transducers:

The COM-DEX is using a micro USB cable for charging.

Guidance and manufacturer's declaration

Electromagnetic emissions

The COM-DEX is intended for use in the electromagnetic environment specified below. The customer or the user of a COM-DEX should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The COM-DEX uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The COM-DEX is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable *)	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable *)	

*) *Battery powered equipment*

Electromagnetic immunity

The COM-DEX is intended for use in the electromagnetic environment specified below. The customer or the user of a COM-DEX should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electro-static discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

<p>Electrical fast transients/ burst IEC 61000-4-4</p>	<p>± 2 kV for power line supplies ± 1 kV for input/output lines</p>	<p>Not applicable *) Not applicable *)</p>	<p>Not applicable *)</p>
<p>Surge IEC 61000-4-5</p>	<p>± 1 kV line(s) to line(s) ± 2 kV line(s) to earth</p>	<p>Not applicable *) Not applicable *)</p>	<p>Not applicable *)</p>

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p><5 % U_T (>95 % dip in U_T) for 0.5 cycle</p> <p>40 % U_T (60 % dip in U_T) for 5 cycles</p> <p>70 % U_T (30 % dip in U_T) for 25 cycles</p> <p><5 % U_T (>95 % dip in U_T) for 5 s</p>	<p>Not applicable *)</p>	<p>Not applicable *)</p>
----------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------	--------------------------

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
---------------------------------------------------------------	-------	-------	------------------------------------------------------------------------------------------------------------------------------------------

NOTE U_T is the a.c. mains voltage prior to the application of the test level.

**) Battery powered equipment*

Electromagnetic immunity – cont.

The COM-DEX is intended for use in the electromagnetic environment specified below. The customer or the user of a COM-DEX should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the COM-DEX, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Conducted RF	3 Vrms	3 Vrms	Recommended separation distance $d = 1.2 \sqrt{P}$
IEC 61000-4-6	150 kHz to 80 MHz		
Radiated RF	3 V/m	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz
IEC 61000-4-3	80 MHz to 2.5 GHz		$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the COM-DEX is used exceeds the applicable RF compliance level above, the COM-DEX should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or re-locating the COM-DEX.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances

Recommended separation distances between portable and mobile RF communication equipment and the COM-DEX.

The COM-DEX is intended for use in the electromagnetic environment in which RF disturbances are controlled. The customer or the user of the COM-DEX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the COM-DEX as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73

1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

This COM-DEX may be interfered with by other equipment even if that other equipment complies with CISPR emission requirements.

EMC/EMI Compliance

The COM-DEX complies with the following EMC/EMI standards:

Standard	Test type	Note
47 CFR Part 15, subpart C	RF emissions	USA Federal Communications Commission (FCC) requirements to intentional radiators.
EN 300 330-2 V1.6.1	RF emissions incl. Spurious emission	EMC and radio spectrum matters for Short Range Devices in the frequency range 9 kHz – 25 MHz
EN 300 328 V1.9.1	RF emissions incl. Spurious emission	EMC and radio spectrum matters for data transmission systems in the 2.4 GHz ISM band and using wide band modulation techniques
IEC 60601-1-2:2007	EMC emission Immunity, RF and ESD	Medical electrical equipment. General requirements for basic safety and essential performance. Electromagnetic compatibility.

EN 301 489-3 V1.6.1	Immunity, RF and ESD	Standard for Low Power Transmitters in the frequency range 9 kHz – 40 GHz
EN 301 489-17 V2.2.1	Immunity, RF and ESD	Standard for Broadband Data Trans- mission Systems in the frequency range 9 kHz – 40 GHz

Hereby, Widex A/S declares that this COM-DEX is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC.

A copy of the Declaration of Conformity can be found at:
<http://www.widex.com/doc>

SYMBOLS

Symbols commonly used by Widex A/S in medical device labeling (labels/IFU/etc.)

Symbol	Title/Description
	Manufacturer The product is produced by the manufacturer whose name and address are stated next to the symbol. If appropriate, the date of manufacture may also be stated.
	Date of manufacture The date when the product was manufactured.
	Use-by date The date after which the product is not to be used.
	Batch code The product's batch code (lot or batch identification).

Symbol	Title/Description
	<p>Catalog number The product's catalog (item) number.</p>
	<p>Serial number The product's serial number.*</p>
	<p>Keep away from sunlight This product must be protected from light source and must be kept away from heat.</p>
	<p>Keep dry The product must be protected from moisture and/or The product must be kept away from rain.</p>
	<p>Lower limit of temperature The lowest temperature to which the product can be safely exposed.</p>

Symbol	Title/Description
	Upper limit of temperature The highest temperature to which the product can be safely exposed.
	Temperature limits The highest and lowest temperatures to which the product can be safely exposed.
	Consult instructions for use The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.
	Caution/Warning Text marked with a caution/warning symbol must be read before using the product.

Symbol	Title/Description
	<p>WEEE mark “Not for general waste”</p> <p>When the product is to be discarded, it must be sent to a designated collection point for recycling and recovery.</p>
	<p>CE mark</p> <p>The product is in conformity with the requirements set out in European CE marking directives.</p>
	<p>Alert</p> <p>The product is identified by R&TTE Directive 1999/5/EC as an equipment Class 2 product with some restrictions on use in some CE member states.</p>
	<p>RCM mark</p> <p>The product complies with EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand market.</p>

Symbol	Title/Description
	Interference Electromagnetic interference may occur in the vicinity of the product.

*The six- or seven-digit number on the product is the serial number. Serial numbers may not always be preceded by **SN**

FCC AND IC STATEMENTS

FCC ID: TTY-CMDEX

IC: 5676B-CMDEX

Federal Communications Commission Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the

equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications to the equipment not expressly approved by Widex could void the user's authority to operate the equipment.

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Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada.

To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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WIDEX A/S